

Coverage Analysis for NGS Data Experiments

Marcel Caraciolo, Bioinformatician and CTO marcel@genomika.com.br




- Introduction to Coverage Analysis
- Main metrics and evaluation methods
- Hands-on (GATK, BedTools, Chanjo)
- Conclusions and further steps



What's Coverage?

Coverage

The read depth in a single position or how many times the base has been uniquely interrogated.






What's Coverage?

As for a set of bases, such as a target region, it represents the average read depth across the given interval.

Average DP: 5.66

DP2 DP3 DP3 DP3 DP3 DP3 DP3







What's Coverage?

As for a set of bases, such as a target region, it represents the average read depth across the given interval.

Average DP: 5.66

DP2 DP3 DP3 DP3 DP3 DP3 DP3







Coverage Analysis

Targeted enrichment sequencing by NGS is a common approach to interrogate specific loci or the whole exome, for instance, in the human exome.

Main problems:

- The efficiency and the lack of bias in the enrichment process needs to be assessed as a quality control step prior to downstream analysis of sequencing data.
- It can be used to distinguish and discard random sequencing and alignment errors. Disease-causing variants can likewise be missed due to poor coverage.
- Save your time, and of course to avoid incorrect conclusion from the analysis.



Coverage Analysis

Illustrated motivation and example:




Illustration of a genomic region in a heterozygote case. The figure shows patient's two alleles of the investigated gene.

The companion of the identified mutation is missed in an area of poor coverage.



Coverage Analysis

Rna-Seq and transcript analysis

We don't know the size/length of the transcriptome.

The best approach is to perform a pseudo coverage (count of alignments spanning a genomic position)


Metrics: RPKM (Reads per KB per million reads)




Some important definitions:

CDS (coding regions of gene)

Portion of a gene's DNA or RNA, composed of exons, that codes for protein.











Some important definitions:

Loci coordinates (genomic coordinates)


Specific location of a gene, DNA Sequence or position on a chromosome.







Some important definitions:




http://genome.ucsc.edu/cgi-bin/hgGateway



Some important definitions:


UCSC coordinates convention.





Some important definitions:

UCSC coordinates convention.



chr<chromsome>:<coordinate>




What files are important in the coverage analysis?

BAM files

It's the binary version of a SAM file. SAM file is a tab-delimited text file that contains the sequence alignment data.







Hands-on: IGV



It requires Java Runtime at your machine.

http://www.broadinstitute.org/software/igv/download

Open at your IGV: http://www.bioinfomgp.org/_media/ngscat/download/example1.bam



What files are important in the coverage analysis?

BED files

It's a tab-delimited text file that defines a feature track (target regions that we are interested to investigate).

```
chr1 213941196 213942363
chr1 213942363 213943530
chr1 213943530 213944697
chr2 158364697 158365864
chr2 158365864 158367031
chr3 127477031 127478198
chr3 127478198 127479365
chr3 127479365 127480532
chr3 127480532 127481699
```



```
What files are important in the coverage
                                                 EGFR_NM_005228_exon_1_chr7_f=
                 chr7
                         55086970
                                     55087058
                         55289978
                                                EGFR_NM_005228_exon_2_chr7_f
                 chr7
                                     55210130
anal
                                                EGFR_NM_005228_exon_3_chr7_f
                 chr7-
                         55218997
                                     55211181
                         55214298
                                     55214433
                                                EGFR_NM_005228_exon_4_chr7_f
                 chr7
                         55218986
                                     55219055
                                                EGFR_NM_005228_exon_5_chr7_f
                 chr7
                         55220238
                                     55220357
                                                 EGFR_NM_005228_exon_6_chr7_f
                         55221703
                                     55221845
                                                EGFR_NM_005228_exon_7_chr7_f
                 chr7
                                                EGFR_NM_005228_exon_8_chr7_f
                 chr7
                         55223522
                                     55223639
                         55224225
                                     55224352
                                                 EGFR_NM_005228_exon_9_chr7_f
                 chr7
                                                EGFR_NM_005228_exon_10_chr7_f
                         55224451
                 chr7
                                     55224525
                 chr7
                         55225355
                                     55225446
                                                EGFR_NM_005228_exon_11_chr7_f
                 chr7-
                         55227831
                                     55228031
                                                 EGFR_NM_005228_exon_12_chr7_f
                 chr7-
                         55229191
                                     55229324
                                                EGFR_NM_005228_exon_13_chr7_f
                 chr7
                         55231425
                                     55231516
                                                 EGFR_NM_005228_exon_14_chr7_f
                 chr7
                         55232972
                                     55233130
                                                 EGFR_NM_005228_exon_15_chr7_f
                         55238867
                                                 EGFR_NM_005228_exon_16_chr7_f
                 chr7
                                     55238906
                         55248675
                 chr7
                                     55240817
                                                EGFR_NM_005228_exon_17_chr7_f
                 chr7
                         55241613
                                     55241736
                                                 EGFR_NM_005228_exon_18_chr7_f+
                         55242414
                                     55242513
                                                 EGFR_NM_005228_exon_19_chr7_f
                 chr7
                 chr7
                         55248985
                                     55249171
                                                EGFR_NM_005228_exon_20_chr7_f
                                                                                                   63
                         55259411
                                                 EGFR_NM_005228_exon_21_chr7_f*
                 chr7
                                     55259567
                         55268458
                                                 EGFR_NM_005228_exon_22_chr7_f
                 chr7
                                     55260534
                                                                                                   30
                 chr7
                         55266409
                                     55266556
                                                 EGFR_NM_005228_exon_23_chr7_f
                                                                                                   97
                 chr7
                         55268008
                                     55268106
                                                 EGFR_NM_005228_exon_24_chr7_f>
                                                                                                   64
                 chr7
                         55268880
                                     55269848
                                                 EGFR_NM_005228_exon_25_chr7_f >
                         55269427
                                     55269475
                                                EGFR_NM_005228_exon_26_chr7_f
                                                                                                   31
                 chr7
                 chr7
                         55270209
                                     55270318
                                                 EGFR_NM_005228_exon_27_chr7_f >
                                                                                                   98
                 chr7
                         55272948
                                     55273310
                                                 EGFR_NM_005228_exon_28_chr7_f
                                                                                                   65
                         115258671
                                                NRAS_PROBE_2_chr1_r=0
                 chr1
                                     115258798
                                                                                                   32
                         115256484
                 chr1
                                     115256570
                                                NRAS_PROBE_3_chr1_r+0+
                 chr1
                         115252220
                                     115252348
                                                NRAS_PROBE_4_1_chr1_r
                                                                                                   99
                                                NRAS_PROBE_4_2_chr1_r
                         115252190
                                     115252239
                 chr1
                         25398252
                                     25398315
                                                KRAS_PROBE_2_chr12_r>
                 chr12
                 chr12
                         25380220
                                     25380317
                                                KRAS_PROBE_3_chr12_r>
                                                KRAS_PROBE_4_1_chr12_r 0
                 chr12
                         25378600
                                     25378721
                         25378548
                                     25378615
                                                KRAS_PROBE_4_2_chr12_r> 0
                 chr12-
```



Hands On Bed file

Â	Genomes	Genome Browser	Tools	Mirrors	Downloads	My Data	
Table Browser							
clade:	Mammal	genome: Hur	man	o ass	embly: Feb.	2009 (GRCh37/h	g19) 🗘
group.	Genes and Ge	ene Prediction Tracks	○ trac	k: RefSeq	Genes	0	
add custom tracks track nubs							
table: refGene							
region: ● genome ○ ENCODE Pilot regions ○ position chr21:33031597-33041570							
Lookup define regions							
identifiers (names/accessions): paste list upload list							
filter: create							
intersection: create							
correlation: create							
	format: BED	- browser extensible	data	٥	Send output to	o ☑ Galaxy →)
GREAT (leave block to keep output in browser)							
output file: (leave blank to keep output in browser) file type returned: plain text gzip compressed							
me type retained. So plaintent Supplemblessed							
get out	put	ry/statistics					

http://genome.ucsc.edu/cgi-bin/hgTables



Hands-on: IGV



Open a bed file with your bam at your IGV.

http://www.broadinstitute.org/software/igv/download

Open at your IGV: http://www.bioinfomgp.org/_media/ngscat/download/example1.bam

http://www.bioinfomgp.org/_media/ngscat/download/seqcap.example1.bed





Main metrics

X coverage or depth coverage

Most common coverage metric. For a simple position it means the read depth. For a set of bases, such as target region, it represents the average read depth across the given interval.

% coverage or breadth of coverage


It representes to what degree a given sequence has been covered at a given read depth. A cutoff for acceptable variant calling can be set at 30x read depth. % coverage means what ratio of sequence that passes our quality cutoff.

Physical coverage

It means the structural coverage for the whole chromosome. It is important for reliably calling structural variations.




Main visualisations






Main visualisations





Main visualisations Coverage distribution






Main visualisations



ERR089819.bam



Coverage (X)



Tools

Most used tools are:

Bedtools, open-source, Python and bash

GATK, open-source, Java and bash

Chanjo, open-source, Python and bash

BioConductor, open-source, R and bash



Bedtools

Fast, flexible toolset for genomic arithmetic.

http://bedtools.readthedocs.org/





GATK

Genome Analysis Toolkit to analyze high-throughput sequencing data.



```
java -Xmx2g -jar GenomeAnalysisTK.jar \
    -R ref.fasta \
    -T DepthOfCoverage \
    -o file_name_base \
    -I input_bams.list
[-geneList refSeq.sorted.txt] \
    [-pt readgroup] \
    [-ct 4 -ct 6 -ct 10] \
    [-L my_capture_genes.interval_list]
```

https://www.broadinstitute.org/gatk/gatkdocs/ org_broadinstitute_gatk_tools_walkers_coverage_DepthOfCoverage.php



GATK

Genome Analysis Toolkit to analyze high-throughput sequencing



1 you supply the salata gument, septimore reage in output an additional cultimary inc that room as ionome.

```
Gene_Name
             Total_Cvg
                                          Sample_1_Total_Cvg
                                                              Sample 1 Avg Cvg
                                                                                  Sample_1_Cvg_Q3
                            Avg_Cvg
                                                                                                       Sample
1_Cvg_Median
                 Sample_1_Cvg_Q1
        594710 238.27 594710 238.27 165
SORT1
                                             245
                                                      330
                                                     >500
NOTCH2 3011542 357.84 3011542 357.84 222
                                             399
LMNA
       563183 186.73 563183 186.73 116
                                             187
                                                     262
NOS1AP 513031 203.50 513031 203.50
                                                     290
                                             191
```

```
[-ct 4 -ct 6 -ct 10] \
[-L my_capture_genes.interval_list]
```

https://www.broadinstitute.org/gatk/gatkdocs/org_broadinstitute_gatk_tools_walkers_coverage_DepthOfCoverage.php



Chanjo

Coverage analysis for clinical sequencing.



```
$ cat intervals.bed | chanjo annotate alignment.bam
#{"sample_id": "bavewira", ...}
1    10    15    interval-1    9.922    0.97231
2    45    55    interval-2    14.231    1.0
```


http://www.chanjo.co/en/latest/index.html



Bioconductor

http://master.bioconductor.org/help/course-materials/2010/SeattleJan10/day3/CoverageEDA.pdf







Conclusions

Simple introduction to Coverage Analysis

Main tools for downstream analysis

Further metrics: GC content, specificity, sensitivity, uniformity.



References

Asan, Xu, Y., Jiang, H., Tyler-Smith, C., Xue, Y., Jiang, T., Wang, J., Wu, M., Liu, X., Tian, G., Wang, J., Wang, J., Yang, H., and Zhang, X. (2011). Comprehensive comparison of three commercial human whole-exome capture platforms. Genome Biology, 12(9), R95.

Clark, M. J., Chen, R., Lam, H. Y. K., Karczewski, K. J., Chen, R., Euskirchen, G., Butte, A. J., and Snyder, M. (2011). Performance comparison of exome DNA sequencing technologies. Nature Biotechnology, 29(10), 908–914.

Frommolt, P., Abdallah, A. T., Altmller, J., Motameny, S. e., Thiele, H., Becker, C., Stemshorn, K., Fischer, M., Freil inger, T., and Nrnberg, P. (2012). Assessing the enrichment performance in targeted resequencing experiments. Human Mutation, 33(4), 635–641.

Hummel, M., Bonnin, S., Lowy, E., and Roma, G. (2011). TEQC: an r package for quality control in target capture experiments. Bioinformatics, 27(9), 1316–1317. Quinlan, A. R. and Hall, I. M. (2010). BEDTools: a flexible suite of utilities for comparing genomic features. Bioinformatics, 26(6), 841–842.

Sikkema-Raddatz, B., Johansson, L. F., de Boer, E. N., Almomani, R., Boven, L. G., van den Berg, M. P., van Spaendonck-Zwarts, K. Y., van Tintelen, J. P., Sijmons, R. H., Jongbloed, J. D. H., and Sinke, R. J. (2013). Targeted next-generation sequencing can replace sanger sequencing in clinical diagnostics. Human Mutation, 34(7), 1035–1042.

Tewhey, R., Nakano, M., Wang, X., Pabon-Pe´na, C., Novak, B., Giuffre, A., Lin, E., Happe, S., Roberts, D. N., LeProust, E. M., et al. (2009). Enrichment of sequencing targets from the human genome by solution hybridization. Genome Biol, 10(10), R116



Coverage Analysis for NGS Data Experiments

Marcel Caraciolo, Bioinformatician and CTO marcel@genomika.com.br