


### Introdução ao NGS


Rodrigo Bertollo rodrigo@genomika.com.br

## Pipeline





# **Terminology**





# Types of Illumina fragment libraries

## single-end




independent reads

### paired-end



two inwardly oriented reads separated by ~200 nt

### mate-paired



two outwardly oriented reads separated by ~3000 nt

# Read Sequences Quality Control

#### Garbage in = garbage out

- Contaminated with other samples?
- Sample barcodes removed?
- Adaptor sequences trimmed?
  - RNAseq, MiSeq data
- Trim ends of reads with poor quality?
  - de novo Assembly
- Know your data
  - Paired reads? Relative orientations?
  - Technology specific concerns?
    - Indels with 454





## Read Sequences

#### **FASTQ Format**

```
@HWI-EAS216_91209:1:2:454:192#0/1
GTTGATGAATTTCTCCAGCGCGAATTTGTGGGCT
+HWI-EAS216_91209:1:2:454:192#0/1
B@BBBBBBBBBBBBAAAA>@AABA?BBBBAAB??>A?
```

Line 1: @read name

Line 2: called base sequence

Line 3: +read name (optional after +)

Line 4: base quality scores



### **FASTQ** format

Standard Format for NGS data

Conversion can be done from sff, fasta + qual, . . .

Extension of the Fasta format

Text-based formats (easy to use!)

If not compressed, it can be huge <a href="http://en.wikipedia.org/wiki/FASTQ\_format">http://en.wikipedia.org/wiki/FASTQ\_format</a>



# Decipher base quality scores

#### Probability of Error = $10^{-Q/10}$

(This is a **Phred** score, also used for other types of qualities.)

- \* Very low quality scores can mean something special Illumina Q ≤ 3 means something like: "I'm lost, you might want to stop believing sequencing cycles from here on out."
- \* In older FASTQ files, the formula and ASCII offset might differ. Consult: <a href="http://en.wikipedia.org/wiki/FASTQ">http://en.wikipedia.org/wiki/FASTQ</a> format



### **FASTQC**

#### **Quality Assurance tool for FASTQ sequences**

FastQC website:

http://www.bioinformatics.babraham.ac.uk

FastQC report documentation:

http://www.bioinformatics.babraham.ac.uk/projects/fastqc/Help/3%20Analysis%20Modules/

Good Illumina dataset:

http://www.bioinformatics.babraham.ac.uk/projects/fastqc/good\_sequence\_short\_fastqc/fastqc\_report.html

Bad Illumina dataset:

http://www.bioinformatics.babraham.ac.uk/projects/fastqc/bad\_sequence\_fastqc/fastqc\_report.html



### **FASTQC**

- Report basic statistics on your data
- Identify issues with your data

#### Basic Statistics

| Measure            | Value                   |  |  |
|--------------------|-------------------------|--|--|
| Filename           | tmp.fastq               |  |  |
| File type          | Conventional base calls |  |  |
| Encoding           | Illumina 1.5            |  |  |
| Total Sequences    | 250000                  |  |  |
| Filtered Sequences | 0                       |  |  |
| Sequence length    | 101                     |  |  |
| %GC                | 51                      |  |  |

#### **<b>№**FastQC Report

#### Summary

- Basic Statistics
- Per base sequence quality
- Per sequence quality scores
- Per base sequence content
- Per base GC content
- Per sequence GC content
- Per base N content
- Sequence Length Distribution
- Sequence Duplication Levels
- Overrepresented sequences
- Kmer Content




## Useful reports

- Should I trim low quality bases?
  - Per-base sequence quality Report
    - based on all sequences
- Do I need to remove adapter sequences?
  - Overrepresented sequences Report
    - based on 1<sup>st</sup> 200,000 sequences
- How complex is my library?
  - Sequence duplication levels Report
    - estimate based on 1<sup>st</sup> 200,000 sequences




# Per Base Sequencing Quality

Quality scores across all bases (Illumina >v1.3 encoding)






# Per Sequencing Quality






# Per Base Sequencing content






### Per Base GC Content




# Per Sequencing Nucleotide Content






### Per Base N content






# Sequence Length Distribution






# Duplicate Sequences Distribution





## Overrepresented K-mers



## Overrepresented K-mers


- What is a k-mer?
- Create a sliding window of size k, move it over all your reads and count occurrence of k-mers
- We can use this to correct sequencing errors!

DNA: ACGTGTAACGTGACGTTGGA

ACGTG


CGTGT

GTGTA





## Overrepresented K-mers




Concept: Rare k-mers are seq. errors Need > 15X coverage


ACGTGGTTACCCTTAAA



# Sequencing Process: PCR primers



# PCR primers





# NGS adaptors & Cutadpt





# NGS adaptors & Cutadpt

Overrepresented sequences

| Sequence   | Count   | Percentage          | Possible Source                          |
|--|---------|---------------------|--|
| GATCGGAAGAGCACACGTCTGAACTCCAGTCACATCACGATATCGTATGC | 1547768 | 38.192098035156306  | TruSeq Adapter, Index 1 (98% over 50bp)  |
| GATCGGAAGAGCACACGTCTGAACTCCAGTCACATCACGATCTCGTATGC | 146635  | 3.61830603513262    | TruSeq Adapter, Index 1 (100% over 50bp) |
| SATCGGAAGAGCACACGTCTGAACTCCAGTCACATCAAGATATCGTATGC | 6639    | 0.16382128255358863 | TruSeq Adapter, Index 1 (97% over 41bp)  |
| SATCGGAAGAGCACACGTCTGAACTCCAGTCACATCACGATTTCGTATGC | 6462    | 0.15945370204267054 | TruSeq Adapter, Index 1 (98% over 50bp)  |
| GATCGGAAGAGCACACGTCTGAACTCCAGTCACATTACGATATCGTATGC | 5433    | 0.1340625136486891  | TruSeq Adapter, Index 1 (97% over 41bp)  |
| ATCGGAAGAGCACACGTCTGAACTCCAGTCACATAACGATATCGTATGC  | 5147    | 0.1270052931621209  | TruSeq Adapter, Index 1 (97% over 41bp)  |
| GATCGGAAGAGCACACGTCTGAACTCCAGTCACACCACGATATCGTATGC | 4703    | 0.11604932849066535 | TruSeq Adapter, Index 1 (97% over 41bp)  |

Very important if your DNA fragment is shorter than read length



# Coverage

 Coverage/depth is how many times that your data covers the genome (on average)

#### Example:

N: Number of reads: 5 mill

L: Read length: 100

G: Genome size: 5 Mbases

$$\bullet$$
 C =  $5*100/5 = 100X$ 

On average there are 100 reads covering each position in the genome

$$C = N \times \frac{L}{G}$$



### Pré-processando os FASTQ's

Marcel Caraciolo, CTO marcel@genomika.com.br